

Digital Video 3D Microscope

JBM- MZ3Di




- An all-in-one digital inspection system for industrial and laboratory use, complete with high-definition imaging, a large 11.6" screen, and a 144-LED ring light.
- Includes a unique 3D rotary attachment which allows you to view object from various angles
- 0.7X-5X zoom magnification is ideal for dissection, soldering, and general inspection, and detents at key magnifications ensure ease and consistency.

High-Accuracy 3D Microscope Images with the Accurate D.F.D. Method



Product Features

Brand Name	JB MICROSCOPE
MODEL	JBM-MZ3Di
Optical System	Finite-conjugate
Objective Magnification	0.7X-5X zoom with detent
Integrated Reduction	0.4X
Working Distance	100mm
Sensor	Sony IMX307 (color)
Sensor Type	CMOS
Sensor Optical Format	1/2.8"
Active Pixels	2M (1920 x 1080)
Active Sensor Area	5.57mm x 3.13mm
Display Size	11.6"
Magnification	Monitor : 46X, Displayed : 12.9X-92X
Native Resolution	1920 x 1080
Illuminator Type	144-LED ring light
Power Requirements	12VDC 1A input, 110-240VAC, 50-60HZ adapter included
Stand Type	Table stand with pillar
Base Dimensions	320mm x 260mm x 20mm





J-9/36

RAJOURIGARDEN

NEW DELHI 110027

www.jbmicroscope.com

info@jbinstruments.com

+91 9354741245

