


Stereo Zoom Microscope With Clamp Table Type JBM-50B



Specification:

Brand	JB Microscope
Optical system	slanting path of rays changes a time of system, stereoscope angle 13°, 45° inclined observation
Objective	0.6X to 50X
Eyepiece	WF10X/20mm Dioptres Adjustment Eyepiece
Base type	table clamp flex arm stand microscope
Magnification	0.6X to 50X
Interpupillary distance adjustment range	55-75mm
Working distance	105mm
Pillar high	240mm
Diameter	φ32mm
Base Size	285x238mm
Type	Optical
Light source	56 LED ring Lights Adjustable brightness
Colour	White+ Black
Product name	Binocular Microscope Stereo Microscope
Objective	0.6X to 50X
Eyepiece	WF10X/20mm





J-9/36

RAJOURIGARDEN

NEW DELHI 110027

www.jbmicroscope.com | info@jbinstruments.com | +91 9354741245