


Stereo Zoom Microscope

JBM-85R4L



Product Specifications

BRAND	JB MICROSCOPE
Model Number	JBM-85R4L
Binocular Head	45 degree inclined 360° rotatable
Interpupillary Distance Adjustment:	55-75mm
Diopter Adjustment	±5dp
Eyepiece	WF10x/20mm
Zoom Objective	0.7x-45x
Zoom Ratio	1:6.4
Magnification	7x-45x
Working Distance	100mm
Head Mounting Size	76mm in diameter
Dust cover	Yes (Included)
Vertical post height	240mm
Vertical post diameter	φ32mm
Built-in wide voltage design	110v~240v
Base size	205x275x40mm
Material	Metal



+++++



J-9/36

RAJOURIGARDEN

NEW DELHI 110027

www.jbmicroscope.com

info@jbinstruments.com

+91 9354741245

+++++