

STEREO ZOOM MICROSCOPE JBM-745B





JBM-745B Series Stereo Zoom Microscope can provide crystal clear sharp images with high resolution and good depth within a broad field of view. Once brought in focus, images remains clear and sharp through out the whole zoom range. It is widely used in electronic industry, educational demonstration and agricultural research for micro-object observation, inspection and measurement.

Brand Name	JB MICROSCOPE
MODEL	JBM-745B
Magnification	45x
Angle	45 Degree
Eyepiece	Standard: 10X/22mm, high eyepoint, wide angle, wide field eyepiece
Objective	Zoom Range 0.67X-4.5X (Standard),
Zoom Ratio	6.7:1, ensure parfocality of the image plane
Illumination	Fluorescent lamp, LED ringlight, for more option.
Aux Lens	0.5X, 0.7X, 0.75X, 1.5X, 2X
Focusing Mount	A1 focusing mount with vertical working distance 50mm
Stand	B1 pole stand, for more options.

[illegible]

+++++



J-9/36

RAJOURIGARDEN

NEW DELHI 110027

www.jbmicroscope.com

info@jbinstruments.com

+91 9354741245

+++++