

Stereo Zoom Microscope

JBM-760T



Brand Name	JB MICROSCOPE
MODEL	JBM-760T
Magnification	7X-60X Stereo Zoom. 3.5X-200X with Optional lens
Eyepiece	WF10X/22mm high eye- Point Wide field eyepiece, convenience for eye -wearing observers WF 15X/15mm, WF20 X10mm, WF25x/9mm eyepiece is optional.
Microscope Head	Binocular Trinocular head, 45°inclined, 360° rotatable. interpupillary distance: 54-76mm. Diopter adjustable on both tubes.
Objective	0.7X-6X stereo zoom, Zoom Ratio: 8.6:1
Auxiliary Objective	0.5X, 0.3X, 0.7X, 2X
Stands	various stand can be configured based on requirement for more options, please refer to page "microscope stand"
illumination:	5V3W up-down led illuminated/Led ring lamp, Refer to page "External illumination"
Camera adapter	phototube, C -mount for Camera

Eyepiece	Standard Configuration Magnification	Working Distance 120mm field of view (mm)
WF10X/22mm	7x 60x	31.4 3.7
WF20X10mm	14X 20X	14.3 2

[illegible]

++++++
++++++
++++++
++++++



J-9/36

RAJOURIGARDEN

NEW DELHI 110027

www.jbmicroscope.com

info@jbinstruments.com

+91 9354741245

++++++
++++++
++++++
++++++