

Binocular Microscope

JBM-SM03





Product Specifications

BRAND	JB MICROSCOPE
Model Name/Number	JBM-SM03
Eye Pieces	10x 18mm
Interpupillary Distance Adjustable	37.7mm
Double Layer Graduated Mechanical Rectangular Stage Size	Double plate mechanical stage size 140mm X 132mm having low positioned co-axial controls
Optional Swing Out Abbe Condenser	yes
Interpupillary distance	48mm to 75mm

- Inter-pupillary distance from 48mm to 75mm. Eyepiece (Anti fungus coated
- S-Plan Achromatic 4X, 10X, 40X, 100X Oil(S.L.) Mechanical Body
- Quadruple nosepiece revolves on ball bearing
- Double plate mechanical stage size 140mm X 132mm
- Cross movement 78mm X 54mm. Sub Stage Condenser
- Abbe condenser N.A.1.25 swing out filter holder, iris diaphragm Illumination System.

Built-in-base illumination system with 3W LED bulb and adjustable brightness.

[illegible]

++++++
++++++
++++++
++++++



J-9/36

RAJOURIGARDEN

NEW DELHI 110027

www.jbmicroscope.com

info@jbinstruments.com

+91 9354741245

++++++
++++++
++++++
++++++