


Binocular Microscope JBM-SM05



Product Specifications

BRAND	JB MICROSCOPE
Model Name/Number	JBM-SM05
Stand	Sturdy stand anti rust materials. Extended base. Viewing Bodies- 30° inclined, 360° rotatable dovetail, Inter papillary distance 54 – 74mm
Eyepiece	Wide Field eyepiece WF10x/18mm, with foldable eye guard, anti-fungus coating – locked in ocular tube.
Nosepiece	Quadruple nosepiece (Ball bearing type)
Objectives	Semi-Plan Achromatic objectives 10x, 40x (spring loaded), 100x Oil (spring loaded) anti-fungus coating. Par-centered & Para focal.
Mechanical Stage	Mechanical stage with Co-axial drive, highly smooth on ball slides, stainless steel holder for slide movement X-Y = 75mm x 50mm.



+++++



J-9/36

RAJOURIGARDEN

NEW DELHI 110027

www.jbmicroscope.com

info@jbinstruments.com

+91 9354741245

+++++