


Biological Microscope

JBM-201B





Technical Parameters

Model Name	JBM-201B
Optical system	Finite Color, infinite Color Corrected Optical Syste
Nosepiece	Quadruple nosepiece
Focus adjustment	Integrated design, Coaxial focus system with upper limited and tension adjustment, Coarse rang :25mm: fine precision :0.002mm: provided with an adjustable tightness device to prevent slipping and a mechanical upper limit device.
Stage	142mm*132mm square, 155mm*142mm round double layers stage, moving range: 76*50mm, Precision: 0.1mm
Condenser	N.A.1.25ABBECCondenserand iris Diaphragm
illumination	3W intensity adjustable LED, bulb center preset, wide voltage 100v-240v,50/60HZ, output DC12V2A
Others	Transmission Filter: blue, yellow, green photo tube: 0.35 x ,0.5x 1x C-mount for camera
viewing Head	30°inclined, gemel binocular head, interpupillary distance: 50-75mm, Diopter ±5 adjustable 30°inclined, gemel trinocular head, interpupillary distance :50-75mm,Diopter ±5 adjustable Splitting Ratio R: T: 80:20
Eyepiece	WF10X/18mm, WF10X/20mm, p110X/18mm optiona
Objective	Achromatic Objective (4x,10x,40x,100x) Plan Achromatic Objective (4x,10x,40x,100x)

++++
++++
++++
++++



J-9/36

RAJOURIGARDEN

NEW DELHI 110027

www.jbmicroscope.com

info@jbinstruments.com

+91 9354741245

++++
++++
++++
++++