

BINOCULAR BIOLOGICAL MICROSCOPE

JBM-B200



- Outstanding image quality, is helpful for pathologist to understand and studied for getting insight of the concerned biological specimen .JBM-B200 Biological Microscope System design for pathological & clinical application,with easy operation and clear image .its ergonomics design, handy coaxial coarse and fine adjustment,meet the demand of pathologist and clinic analyst. micro world, new vision. with PL10X/18MM wide field plan eyepiece,JBM-B200 Creates more comfortable observation experience. combined with optional plan achromatic objectives,JBM-B200 is able to provide a sharp image in natural optional plan achromatic objectives,JBM-B200 is able to provide a sharp image in natural color,immersive yourself to the micro world and satisfy your desire to explore the unknown.

High-Accuracy 3D Microscope Images with the Accurate D.F.D. Method



Technical Parameters

Brand	JB Microscope
Model no	JBM-B200
VIEWING HEAD	30 DEGREE INCLINED BINOCULAR HEAD, WITH LOCKING MECHANISM FOR EYEPIECES, RIGHT DIOPTER
ADJUSTMENT	:±; INTERPUPILLARY DISTANCE RANGE:50-75MM.
EYEPIECE	WIDE FIELD PLAN EYEPIECE PL10X/18MM, IS ABLE TO ATTACH POINTER OR RETICLE.
OBJECTIVE	PLAN ACHROMATIC OBOECTIVES 4X/10X/40X/100X (OIL
FRAME	LOW-POSITION COARSE AND FINE ADJUSTMENT WITH STOP AND TIGHTEN STRUCTURE, COARSE RANGE:13MM (COMBINE WITH FIXED STAGE), 7MM (COMBINE WITH MECHANICAL STAGE);FINE PRECISION 0.002MM; 0.2W LED ILLUMINATION (COLOR TEMPERATURE 5700-6500K), INTENSITY ADJUSTABLE; POWER OUTPUT 100V~240V; POWER INPUT 5V1A; ADJUSTABLE; POWER OUTPUT 100V~240V; POWER INPUT 5V1A; (RECHARGEABLE FRAME FOR OPTION)
NOSEPIECE	: QUADRUPLE NOSEPIECE
STAGE	FIXED STAGE WITH TWO CLIPS, 110MMX120MM; N.A.0.65 SIMPLE CONDENSER WITH FIVE-HOLE ROTATABLE DIAPHRAGM. DOUBLE-LAYER MECHANICAL STAGE, 155MMX125MM; N.A.1.25 ABBE CONDENSER WITH IRIS DIAPHRAGM AND FILTER HILDE

+++++



J-9/36

RAJOURIGARDEN

NEW DELHI 110027

www.jbmicroscope.com

info@jbinstruments.com

+91 9354741245

+++++