


Upright Inspection Microscopes JBM-012R





Product Features

Brand Name	JB MICROSCOPE
MODEL	JBM-012R
Eyepieces	• High Eye-point Wide Field Plan Eyepiece PL10X/25mm, diopter adjustable • High Eye-point Wide Field Plan Eyepiece PL10X/25mm, diopter adjustable, with reticle
Viewing Head	0-35° Tilting trinocular head, Erect image, interpupillary distance: 50-76mm, splitting ratio: 100 : 0 or 0 : 100
Objective	B&D semi-achromatic 5X, 10X, 20X, 50X, 100X
DIC	DIC Slide
Nosepiece	Electronic BD sextuple nosepiece with DIC slot Bright & Dark field quintuple nosepiece (with DIC slot) Bright & Dark field sextuple nosepiece (with DIC slot) Bright field sextuple nosepiece (with DIC slot) Bright field Septuple nosepiece (with DIC slot
Focus	Coaxial Coarse and fine focusing adjustment, Coarse adjustment distance: 35mm, precision of fine adjustment: 0.001mm, with loosen and tight limit ring
Stage	14 X 12 inches three-layer mechanical stage with low position coaxial adjustment, size: 718mm X 420mm, moving range: 356mm X 305mm; with clutch handle for quick movement; glass plate for reflected use.
Condensor	B&D reflect illumination, with Electronic Iris field diaphragm and aperture diaphragm, central adjustable, with filter slot and polarizing slot, with switch device with bright and dark field
Illumination	12V 100W Halogen lamp (for reflection / Transmitted) 12V 100W Halogen lamp (PHILIPS 7724)
Accessories	12 inch Rotatable work plate, apply to 6, 8, 12 inch Wafer.

[illegible]

+++++



J-9/36

RAJOURIGARDEN

NEW DELHI 110027

www.jbmicroscope.com

info@jbinstruments.com

+91 9354741245

+++++