

JBM-SM04



Product Specifications

BRAND	JB MICROSCOPE
Model Number	JBM-SM04
Observation Head	Binocular tube. inclined at 45° rotatable through 360°
Magnification	100x -1000x
Eyepieces	10x
Objectives	5x, 10x, 40x, 100x
Stage	Casting body stage with mechanical stage for fine movement of slides
Focusing	Ergonomic low position co-axial coarse and fine focusing system on ball bearing guide ways.
Condenser	Movable ABBE condenser, with aspheric lens NA 1.25 & daylight blue filter
Illumination	Illumination system 6V-20W Halogen with intensity control. Power Input 220V AC Optional – LED light source

[illegible]

+++++



J-9/36

RAJOURIGARDEN

NEW DELHI 110027

www.jbmicroscope.com

info@jbinstruments.com

+91 9354741245

+++++