


JBM-SME-02





All the Microscopes are latest in design, precision-engineered with finest features, usually found in high-priced models. Their operational functions can be carried out without any difficulty and provide excellent crisp clear image so that even advanced studies are routinely possible.

At high power, you'll be able to examine specimen slides of yeasts and molds, cultures, plant and animal parts, fibers, bacteria and more. Use the lower magnifications to examine objects like coins, rocks, insects and more.



Brand	JB MICROSCOPE
Model Name	JBM-SME-02
Observation Head	Monocular die-cast body, inclinable up to a horizontal position (90°).
Magnification	100x – 1000x
Eyepieces	10X
Objectives	10x, 40x, 60x, 100x
Stage	Casting body stage with clips holding the slides
Focusing	Quadruple Revolving Nosepiece is provided with positive click stop for proper centering
Condenser	Movable ABBE condenser, with aspheric lens NA 1.25 & daylight blue filter
Illumination	Illumination system 6V-20W Halogen with intensity control. Power Input 220V AC

+++++



J-9/36

RAJOURIGARDEN

NEW DELHI 110027

www.jbmicroscope.com

info@jbinstruments.com

+91 9354741245

+++++