

DIGITAL MICROSCOPE CAMERA

JBM-500
(USB 3.0) 5MP



5MP
High Resolution



USB 3.0
High Speed



Ultra HD
Image Quality



Wide
Compatibility



USB2.0 high speed transmission interface
simple and convenient at the same time to
ensure stable high transmission rate full
resolution can realize 30 frames fast
imaging



Technical Specification

Brand	JB MICROSCOPE
Model	JBM-500 (USB 3.0) 5MP
Camera pixel	5 Mega Pixel
Chip mode	IMX335
Optical Size	1/2.8 inch
Image size	2.0 μ mX2.0 μ m
Resolution	2952(H)X 1944(V)
Camera frame rate	up to 30 fps
Output Interface	USB2.0high speed output
External Device	Desktop PC or PC Notebook
Optical Interface	Standard C connector
Auto Settings	Exposure, Brightness, White Balance
Manual Settings	Exposure, Gain, Noise Reduction, Gamma



J-9/36

RAJOURIGARDEN

NEW DELHI 110027

www.jbmicroscope.com

info@jbinstruments.com

+91 9354741245